

Product Datasheet

Human Phospholipase A2 ELISA Kit (orb1214389)

Description	Human Phospholipase A2, PL-A2 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	PLA2G1B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,360
Uniprot ID	P04054
Solubility (25°C)	64ng/ml
Sensitivity	0.096ng/ml
Entrez	5319
Expiration Date	Please enquire.