

Product Datasheet

Human Matrix Metalloproteinase 10 ELISA Kit (orb1214394)

Description	Human Matrix Metalloproteinase 10, MMP-10 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	MMP10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,151
Uniprot ID	P09238
Solubility (25°C)	48ng/ml
Sensitivity	0.06ng/ml
Entrez	4319
Expiration Date	Please enquire.