

Product Datasheet

Human Heme Oxygenase 2 ELISA Kit (orb1214398)

Description	Human Heme Oxygenase 2, HO-2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	HMOX2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,033
Uniprot ID	P30519
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	3163
Expiration Date	Please enquire.