

Product Datasheet

Human Kallikrein 8 ELISA Kit (orb1214405)

Description	Human Kallikrein 8, KLK8 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	KLK8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,048
Uniprot ID	O60259
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	11202
Expiration Date	Please enquire.