

Product Datasheet

Human B-galactosidase ELISA Kit (orb1214406)

Description	Human B-galactosidase, BGAL ELISA Kit
Reactivity	Human
Range	3-900uIU/mL
Tested Applications	ELISA
Target	GLB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	76,075
Uniprot ID	P16278
Solubility (25°C)	960µIU/mL
Sensitivity	1.26µIU/mL
Entrez	2720
Expiration Date	Please enquire.