

Product Datasheet

Human B-glucosidase ELISA Kit (orb1214407)

Description	Human B-glucosidase, B-GLU ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	GLB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	76,075
Uniprot ID	P16278
Solubility (25°C)	6400ng/L
Sensitivity	10.14ng/L
Entrez	2720
Expiration Date	Please enquire.