

Product Datasheet

Human Sialidase-1 ELISA Kit (orb1214413)

Description	Human Sialidase-1, NEU1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	NEU1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,467
Uniprot ID	Q99519
Solubility (25°C)	48ng/mL
Sensitivity	0.045ng/mL
Entrez	4758
Expiration Date	Please enquire.