

Product Datasheet

Human Cystathionine-gamma-lyase ELISA Kit (orb1214431)

Description	Human Cystathionine-gamma-lyase, CSE ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	CTH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,508
Uniprot ID	P32929
Solubility (25°C)	80ng/ml
Sensitivity	0.12ng/ml
Entrez	1491
Expiration Date	Please enquire.