

Product Datasheet

Human Cyclophilin A ELISA Kit (orb1214443)

Description	Human Cyclophilin A, CYPA ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	PPIA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,012
Uniprot ID	P62937
Solubility (25°C)	240ng/ml
Sensitivity	0.27ng/ml
Entrez	5478
Expiration Date	Please enquire.