

Product Datasheet

Human Suppressors of Cytokine Signaling 3 ELISA Kit (orb1214448)

Description	Human Suppressors of Cytokine Signaling 3, SOCS-3 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	SOCS3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,770
Uniprot ID	O14543
Solubility (25°C)	6400ng/L
Sensitivity	9.02ng/L
Entrez	9021
Expiration Date	Please enquire.