

Product Datasheet

Human Forkhead Box O3 ELISA Kit (orb1214456)

Description	Human Forkhead Box O3, FOXO3 ELISA Kit
Reactivity	Human
Range	0.5-200ng/L
Tested Applications	ELISA
Target	FOXO3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,277
Uniprot ID	043524
Solubility (25°C)	240ng/L
Sensitivity	0.33ng/L
Entrez	2309
Expiration Date	Please enquire.