

Product Datasheet

Human Beta-Arrestin 2 ELISA Kit (orb1214471)

Description	Human Beta-Arrestin 2, ARRB2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ARRB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,106
Uniprot ID	P32121
Solubility (25°C)	3200ng/L
Sensitivity	7.39ng/L
Entrez	409
Expiration Date	Please enquire.