

Product Datasheet

Human Beta-Arrestin 1 ELISA Kit (orb1214472)

Description	Human Beta-Arrestin 1, ARRB1 ELISA Kit
Reactivity	Human
Range	0.1-38ng/ml
Tested Applications	ELISA
Target	ARRB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,066
Uniprot ID	P49407
Solubility (25°C)	40ng/ml
Sensitivity	0.05ng/ml
Entrez	408
Expiration Date	Please enquire.