

Product Datasheet

Human Protein Phosphatase 1 Regulatory Subunit 1A ELISA Kit (orb1214474)

Description	Human Protein Phosphatase 1 Regulatory Subunit 1A, PPP1R1A ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	PPP1R1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,011
Uniprot ID	Q13522
Solubility (25°C)	3200ng/L
Sensitivity	7.17ng/L
Entrez	5502
Expiration Date	Please enquire.