

Product Datasheet

Human Complement Factor B ELISA Kit (orb1214478)

Description	Human Complement Factor B, CFB ELISA Kit
Reactivity	Human
Range	8-2000ug/ml
Tested Applications	ELISA
Target	CFB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	85,533
Uniprot ID	P00751
Solubility (25°C)	2400µg/ml
Sensitivity	4.19µg/ml
Entrez	629
Expiration Date	Please enquire.