

Product Datasheet

Human Complement 3 ELISA Kit (orb1214481)

Description	Human Complement 3, C3 ELISA Kit
Reactivity	Human
Range	0.05-30mg/ml
Tested Applications	ELISA
Target	C3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	187,148
Uniprot ID	P01024
Solubility (25°C)	32mg/ml
Sensitivity	0.025mg/ml
Entrez	718
Expiration Date	Please enquire.