

Product Datasheet

Human Toll-like Receptor 1 ELISA Kit (orb1214483)

Description	Human Toll-like Receptor 1, TLR1 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	TLR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,291
Uniprot ID	Q15399
Solubility (25°C)	80ng/ml
Sensitivity	0.09ng/ml
Entrez	7096
Expiration Date	Please enquire.