

Product Datasheet

Human Toll-like Receptor 2 ELISA Kit (orb1214485)

Description	Human Toll-like Receptor 2, TLR2 ELISA Kit
Reactivity	Human
Range	0.1-40.ng/ml
Tested Applications	ELISA
Target	TLR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	89,838
Uniprot ID	O60603
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	7097
Expiration Date	Please enquire.