

Product Datasheet

Human Toll-like Receptor 9 ELISA Kit (orb1214486)

Description	Human Toll-like Receptor 9, TLR-9 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	TLR9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	115,860
Uniprot ID	Q9NR96
Solubility (25°C)	1600ng/L
Sensitivity	2.53ng/L
Entrez	54106
Expiration Date	Please enquire.