

Product Datasheet

Human Toll-like Receptor 7 ELISA Kit (orb1214493)

Description	Human Toll-like Receptor 7, TLR7 ELISA Kit
Reactivity	Human
Range	10-3500ng/L
Tested Applications	ELISA
Target	TLR7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	120,922
Uniprot ID	Q9NYK1
Solubility (25°C)	4000ng/L
Sensitivity	5.32ng/L
Entrez	51284
Expiration Date	Please enquire.