

Product Datasheet

Human Membrane Co Factor Protein ELISA Kit (orb1214505)

Description	Human Membrane Co Factor Protein, MCP ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	CD46
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,747
Uniprot ID	P15529
Solubility (25°C)	480ng/ml
Sensitivity	0.51ng/ml
Entrez	4179
Expiration Date	Please enquire.