

Product Datasheet

Human Intercellular Adhesion Molecule 3 ELISA Kit (orb1214507)

Description	Human Intercellular Adhesion Molecule 3, ICAM-3 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	ICAM3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,541
Uniprot ID	P32942
Solubility (25°C)	480ng/ml
Sensitivity	0.41ng/ml
Entrez	3385
Expiration Date	Please enquire.