

Product Datasheet

Human Cd163 Antigen ELISA Kit (orb1214512)

Description	Human Cd163 Antigen, CD163 ELISA Kit
Reactivity	Human
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	CD63
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,637
Uniprot ID	P08962
Solubility (25°C)	32ng/ml
Sensitivity	0.052ng/ml
Entrez	967
Expiration Date	Please enquire.