

Product Datasheet

Human C-C Chemokine Receptor Type 4 ELISA Kit (orb1214517)

Description	Human C-C Chemokine Receptor Type 4, CCR4 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	CCR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,403
Uniprot ID	P51679
Solubility (25°C)	12.8ng/ml
Sensitivity	0.021ng/ml
Entrez	1233
Expiration Date	Please enquire.