

Product Datasheet

Human Cluster of Differentiation 3 ELISA Kit (orb1214518)

Description	Human Cluster of Differentiation 3, CD3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	EID3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,168
Uniprot ID	Q8N140
Solubility (25°C)	24ng/ml
Sensitivity	0.01ng/ml
Entrez	493861
Expiration Date	Please enquire.