

Product Datasheet

Human Cluster of Differentiation ELISA Kit (orb1214519)

Description	Human Cluster of Differentiation, CD34 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	KDF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,642
Uniprot ID	Q8NAX2
Solubility (25°C)	16ng/ml
Sensitivity	0.012ng/ml
Entrez	126695
Expiration Date	Please enquire.