

## Product Datasheet

### Human E-Cadherin ELISA Kit (orb1214523)

|                                |                                           |
|--------------------------------|-------------------------------------------|
| <b>Description</b>             | Human E-Cadherin, E-CAD<br>ELISA Kit      |
| <b>Reactivity</b>              | Human                                     |
| <b>Range</b>                   | 1-300ng/ml                                |
| <b>Tested<br/>Applications</b> | ELISA                                     |
| <b>Target</b>                  | CDH1                                      |
| <b>Storage</b>                 | 2-8°C                                     |
| <b>Note</b>                    | For research use only                     |
| <b>Application notes</b>       | Sandwich                                  |
| <b>Sample Types</b>            | Serum, plasma, cell culture<br>supernates |
| <b>MW</b>                      | 97,456                                    |
| <b>Uniprot ID</b>              | <b>P12830</b>                             |
| <b>Solubility (25°C)</b>       | 400ng/ml                                  |
| <b>Sensitivity</b>             | 0.52ng/ml                                 |
| <b>Entrez</b>                  | <b>999</b>                                |
| <b>Expiration Date</b>         | Please enquire.                           |