

Product Datasheet

Human Lymphotoxin-beta ELISA Kit (orb1214526)

Description	Human Lymphotoxin-beta, LTB ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	TNFC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,390
Uniprot ID	Q06643
Solubility (25°C)	4800ng/L
Sensitivity	9.98ng/L
Entrez	4050
Expiration Date	Please enquire.