

Product Datasheet

Human Intelectin-1 ELISA Kit (orb1214530)

Description	Human Intelectin-1, ITLN1 ELISA Kit
Reactivity	Human
Range	2-600ng/L
Tested Applications	ELISA
Target	ITLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,962
Uniprot ID	Q8WWA0
Solubility (25°C)	640ng/L
Sensitivity	1.03ng/L
Entrez	55600
Expiration Date	Please enquire.