

Product Datasheet

Human Interleukin 18 ELISA Kit (orb1214533)

Description	Human Interleukin 18, IL-18 ELISA Kit
Reactivity	Human
Range	0.5-100ng/L
Tested Applications	ELISA
Target	IL18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,326
Uniprot ID	Q14116
Solubility (25°C)	128ng/L
Sensitivity	0.2ng/L
Entrez	3606
Expiration Date	Please enquire.