

Product Datasheet

Human Vascular Endothelial Cell Growth Factor B ELISA Kit (orb1214544)

Description	Human Vascular Endothelial Cell Growth Factor B, VEGF-B ELISA Kit
Reactivity	Human
Range	20-1200ng/L
Tested Applications	ELISA
Target	VEGFB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,602
Uniprot ID	P49765
Solubility (25°C)	1280ng/L
Sensitivity	10.11ng/L
Entrez	7423
Expiration Date	Please enquire.