

Product Datasheet

Human Podoplanin ELISA Kit (orb1214548)

Description	Human Podoplanin, PDPN ELISA Kit
Reactivity	Human
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	PDPN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,698
Uniprot ID	Q86YL7
Solubility (25°C)	32ng/ml
Sensitivity	0.05ng/ml
Entrez	10630
Expiration Date	Please enquire.