

Product Datasheet

Human Interleukin 15 ELISA Kit (orb1214559)

Description	Human Interleukin 15, IL-15 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	IL15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,086
Uniprot ID	P40933
Solubility (25°C)	1600ng/L
Sensitivity	1.98ng/L
Entrez	3600
Expiration Date	Please enquire.