

Product Datasheet

Human Interleukin 16 ELISA Kit (orb1214560)

Description	Human Interleukin 16, IL-16 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	IL16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	141752
Uniprot ID	Q14005
Solubility (25°C)	960ng/L
Sensitivity	1.47ng/L
Entrez	3603
Expiration Date	Please enquire.