

Product Datasheet

Human Vascular Endothelial Cell Growth Factor C ELISA Kit (orb1214573)

Description	Human Vascular Endothelial Cell Growth Factor C, VEGF-C ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	VEGFC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,883
Uniprot ID	P49767
Solubility (25°C)	6400ng/L
Sensitivity	1.18ng/L
Entrez	7424
Expiration Date	Please enquire.