

Product Datasheet

Human Vascular Endothelial Cell Growth Factor D ELISA Kit (orb1214574)

Description	Human Vascular Endothelial Cell Growth Factor D, VEGF-D ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	VEGFD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,444
Uniprot ID	O43915
Solubility (25°C)	1280ng/L
Sensitivity	2.51ng/L
Entrez	2277
Expiration Date	Please enquire.