

## Product Datasheet

### Human Interleukin 23 ELISA Kit (orb1214575)

|                            |                                          |
|----------------------------|------------------------------------------|
| <b>Description</b>         | Human Interleukin 23, IL-23<br>ELISA Kit |
| <b>Reactivity</b>          | Human                                    |
| <b>Range</b>               | 3-900ng/L                                |
| <b>Tested Applications</b> | ELISA                                    |
| <b>Target</b>              | IL27                                     |
| <b>Storage</b>             | 2-8°C                                    |
| <b>Note</b>                | For research use only                    |
| <b>Application notes</b>   | Sandwich                                 |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates   |
| <b>MW</b>                  | 27,493                                   |
| <b>Uniprot ID</b>          | <b>Q8NEV9</b>                            |
| <b>Solubility (25°C)</b>   | 960ng/L                                  |
| <b>Sensitivity</b>         | 1.52ng/L                                 |
| <b>Entrez</b>              | <b>246778</b>                            |
| <b>Expiration Date</b>     | Please enquire.                          |