

Product Datasheet

Human Interleukin 25 ELISA Kit (orb1214583)

Description	Human Interleukin 25, IL-25 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	IL25
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,330
Uniprot ID	Q9H293
Solubility (25°C)	6400ng/L
Sensitivity	10.5ng/L
Entrez	64806
Expiration Date	Please enquire.