

## Product Datasheet

### Human Interleukin 26 ELISA Kit (orb1214584)

|                            |                                          |
|----------------------------|------------------------------------------|
| <b>Description</b>         | Human Interleukin 26, IL-26<br>ELISA Kit |
| <b>Reactivity</b>          | Human                                    |
| <b>Range</b>               | 0.5-60ng/L                               |
| <b>Tested Applications</b> | ELISA                                    |
| <b>Target</b>              | IL26                                     |
| <b>Storage</b>             | 2-8°C                                    |
| <b>Note</b>                | For research use only                    |
| <b>Application notes</b>   | Sandwich                                 |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates   |
| <b>MW</b>                  | 19,843                                   |
| <b>Uniprot ID</b>          | <b>Q9NPH9</b>                            |
| <b>Solubility (25°C)</b>   | 64ng/L                                   |
| <b>Sensitivity</b>         | 0.23ng/L                                 |
| <b>Entrez</b>              | <b>55801</b>                             |
| <b>Expiration Date</b>     | Please enquire.                          |