

Product Datasheet

Human Growth Differentiation Factor 2 ELISA Kit (orb1214585)

Description	Human Growth Differentiation Factor 2, GDF2 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	BMP9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,320
Uniprot ID	Q9UK05
Solubility (25°C)	960ng/L
Sensitivity	1.39ng/L
Entrez	2658
Expiration Date	Please enquire.