

Product Datasheet

Human Interleukin 29 ELISA Kit (orb1214593)

Description	Human Interleukin 29, IL-29 ELISA Kit
Reactivity	Human
Range	6-1800ng/L
Tested Applications	ELISA
Target	IL29
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,898
Uniprot ID	Q8IU54
Solubility (25°C)	1920ng/L
Sensitivity	3.42ng/L
Entrez	282618
Expiration Date	Please enquire.