

Product Datasheet

Human Parathyroid Hormone-like Protein ELISA Kit (orb1214595)

Description	Human Parathyroid Hormone-like Protein, PLP ELISA Kit
Reactivity	Human
Range	0.05-10pmol/L
Tested Applications	ELISA
Target	PTH LH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,194
Uniprot ID	P12272
Solubility (25°C)	12.8pmol/L
Sensitivity	0.022pmol/L
Entrez	5744
Expiration Date	Please enquire.