

Product Datasheet

Human Glucagon-like Peptide 1 ELISA Kit (orb1214599)

Description	Human Glucagon-like Peptide 1, GLP1 ELISA Kit
Reactivity	Human
Range	2-180pmol/L
Tested Applications	ELISA
Target	GCG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20909
Uniprot ID	P01275
Solubility (25°C)	192pmol/L
Sensitivity	1.27pmol/L
Entrez	2641
Expiration Date	Please enquire.