

Product Datasheet

Horse Somatotropin ELISA Kit (orb1214605)

Description	Horse Somatotropin, GH1 ELISA Kit
Reactivity	Equine
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	GH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,423
Uniprot ID	P01245
Solubility (25°C)	16ng/ml
Sensitivity	0.017ng/ml
Entrez	100034180
Expiration Date	Please enquire.