

Product Datasheet

Goat Interleukin 18 ELISA Kit (orb1214627)

Description	Goat Interleukin 18, IL-18 ELISA Kit
Reactivity	Goat
Range	10-2000ng/L
Tested Applications	ELISA
Target	IL18
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,023
Uniprot ID	Q3ZT29
Solubility (25°C)	2400ng/L
Sensitivity	5.51ng/L
Entrez	100861190
Expiration Date	Please enquire.