

Product Datasheet

Goat Interleukin 6 ELISA Kit (orb1214628)

Description	Goat Interleukin 6, IL-6 ELISA Kit
Reactivity	Goat
Range	1-300ng/L
Tested Applications	ELISA
Target	IL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,423
Uniprot ID	Q28319
Solubility (25°C)	320ng/L
Sensitivity	0.52ng/L
Entrez	100860785
Expiration Date	Please enquire.