

Product Datasheet

Canine Interleukin-5 ELISA Kit (orb1214658)

Description	Canine Interleukin-5, IL-5 ELISA Kit
Reactivity	Canine
Range	0.5-300ng/L
Tested Applications	ELISA
Target	IL5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,307
Uniprot ID	Q95J76
Solubility (25°C)	320ng/L
Sensitivity	0.28ng/L
Entrez	403790
Expiration Date	Please enquire.