

Product Datasheet

Canine Appetite-Regulating Hormone ELISA Kit (orb1214669)

Description	Canine Appetite-Regulating Hormone, GHRL ELISA Kit
Reactivity	Canine
Range	20-6000ng/L
Tested Applications	ELISA
Target	GHRL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	13,007
Uniprot ID	Q9BEF8
Solubility (25°C)	6400ng/L
Sensitivity	10.55ng/L
Entrez	403587
Expiration Date	Please enquire.