

Product Datasheet

Canine Calcitonin ELISA Kit (orb1214679)

Description	Canine Calcitonin, CT ELISA Kit
Reactivity	Canine
Range	5-1000ng/L
Tested Applications	ELISA
Target	CALCA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,051
Uniprot ID	P41547
Solubility (25°C)	1200ng/L
Sensitivity	2.27ng/L
Entrez	403946
Expiration Date	Please enquire.