

Product Datasheet

Canine Leptin ELISA Kit (orb1214686)

Description	Canine Leptin, LEP ELISA Kit
Reactivity	Canine
Range	2-150ng/ml
Tested Applications	ELISA
Target	LEPTIN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,655
Uniprot ID	002720
Solubility (25°C)	160ng/ml
Sensitivity	1.01ng/ml
Entrez	403616
Expiration Date	Please enquire.